

Anti-MUC2 antibody [EBS-T-233] (STJ16103073)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/IHC-F/IF/WB/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Purification	Affinity purified
Formulation	PBS with 0.02% Sodium Azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	4583
Gene Symbol	MUC2
UniProt ID	MUC2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---