

## Anti-Myeloid Specific Antigen antibody [BM-2] (STJ16103092)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F/IHC-P
Host / Source	Mouse
Reactivity	Human/Macaque Monkey

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---