

Anti-blood group antigen ABH antibody [HE-10] (STJ16103094)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-P FC 1-3:100
Formulation	Culture supernatant with <0.1% sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---