

## Anti-Biotin antibody [Hyb-8] (STJ16103103)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IHC-F/IHC-P/IF/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Rabbit/Species Restrictions Not Known

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---