

Anti-Actin antibody [HHF35] (STJ16103106)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IHC-F/ICC
Host / Source	Mouse
Reactivity	Human/Chicken/Mouse/Non-Human Primates/Dog/Rat/Rabbit/Cat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---