

## Anti-beta-Catenin antibody [EM-22] (STJ16103114)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/WB/ICC
Host / Source	Mouse
Reactivity	Human/Mouse/Hamster

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">1499</a>
Gene Symbol	<a href="#">CTNNB1</a>
UniProt ID	<a href="#">CTNB1_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---