

Anti-beta-Tubulin antibody [TU-12] (STJ16103116)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/FC/WB
Host / Source	Mouse
Reactivity	Mouse/Rat/Pig/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Formulation	Tris buffered saline (TBS) , pH 8.0, 15 mM sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---