

Anti-Clathrin Heavy Chain antibody [BF-06] (STJ16103126)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/WB/ICC/ELISA
Host / Source	Mouse
Reactivity	Rat/Cow/Mouse/Human/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Formulation	Stabilizing Tris buffered saline (TBS) , pH 8.0, 15 mM sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---