

## Anti-HLA-C antibody [DT-9] (STJ16103147)

### GENERAL INFORMATION

|               |                          |
|---------------|--------------------------|
| Product Type  | Primary antibodies       |
| Applications  | IP/FC                    |
| Host / Source | Mouse                    |
| Reactivity    | Non-Human Primates/Human |

### PRODUCT PROPERTIES

|                     |                                                              |
|---------------------|--------------------------------------------------------------|
| Clonality           | Monoclonal                                                   |
| Concentration       | 1 mg/mL                                                      |
| Conjugation         | Unconjugated                                                 |
| Purification        | Purified by protein-A affinity chromatography.               |
| Formulation         | Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide |
| Isotype             | IgG2b                                                        |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.               |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">3107</a>       |
| Gene Symbol | <a href="#">HLA-C</a>      |
| UniProt ID  | <a href="#">HLAC_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|