

## Anti-HLA-DR antibody [MEM-12] (STJ16103148)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC/IP/WB           |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 1 mg/mL  |
| Conjugation         | Unconjugated   |
| Purification        | Purified by protein-A affinity chromatography.               |
| Formulation         | Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide |
| Isotype             | IgG1   |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.               |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|