

Anti-C3aR antibody [HC3aRZ8] (STJ16103157)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC/IHC-F |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Purified by protein-A affinity chromatography. |
| Formulation | Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide |
| Isotype | IgG2b |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 719 |
| Gene Symbol | C3AR1 |
| UniProt ID | C3AR_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|