

Anti-CD170 antibody [1A5] (STJ16103209)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/IHC-F/ICC/ELISA
Host / Source	Mouse
Reactivity	Non-Human Primates/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---