

Anti-CD180 antibody [G28-8] (STJ16103217)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	4064
Gene Symbol	CD180
UniProt ID	CD180_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---