

## Anti-CD184 antibody [12G5] (STJ16103218)

### GENERAL INFORMATION

|               |                          |
|---------------|--------------------------|
| Product Type  | Primary antibodies       |
| Applications  | FC/IHC-P/ICC             |
| Host / Source | Mouse                    |
| Reactivity    | Human/Non-Human Primates |

### PRODUCT PROPERTIES

|                     |                                                              |
|---------------------|--------------------------------------------------------------|
| Clonality           | Monoclonal                                                   |
| Concentration       | 1 mg/mL                                                      |
| Conjugation         | Unconjugated                                                 |
| Purification        | Purified by protein-A affinity chromatography.               |
| Formulation         | Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide |
| Isotype             | IgG2ak                                                       |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.               |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|