

## Anti-CD47 antibody [MEM-122] (STJ16103290)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/WB/IHC-F
Host / Source	Mouse
Reactivity	Non-Human Primates/Human/Pig

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Formulation	Tris buffered saline (TBS) , pH 8.0, 15 mM sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">961</a>
Gene Symbol	<a href="#">CD47</a>
UniProt ID	<a href="#">CD47_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---