

Anti-CD49d antibody [9F10] (STJ16103294)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F
Host / Source	Mouse
Reactivity	Sheep/Cat/Horse/Dog/Cow/Non-Human Primates/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---