

## Anti-CD79a antibody [HM57] (STJ16103318)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-P/IHC-F
Host / Source	Mouse
Reactivity	Chicken/Rabbit/Opossum/Guinea Pig/Horse/Cow/Pig/Mouse/Rat/Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">973</a>
Gene Symbol	<a href="#">CD79A</a>
UniProt ID	<a href="#">CD79A_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---