

Anti-TROP2 antibody [TrMab-6] (STJ16103363)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | WB/FC/IHC/ICC |
| Host / Source | Mouse |
| Reactivity | Pig/Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Purified by protein-A affinity chromatography. |
| Formulation | Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide |
| Isotype | IgG2bk |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|