

Anti-Neurofilament H antibody [NF-05] (STJ16103376)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/ELISA
Host / Source	Mouse
Reactivity	Human/Rat/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Formulation	Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---