

Anti-PCNA antibody [PC10] (STJ16103384)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/WB/IHC-P/IHC-F/ICC
Host / Source	Mouse
Reactivity	Mouse/Human/Non-Human Primates/Drosophila/Chicken/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	5111
Gene Symbol	PCNA
UniProt ID	PCNA_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---