

Anti-RPSA antibody [RP-01] (STJ16103392)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/ELISA/IP/FC
Host / Source	Mouse
Reactivity	Mouse/Xenopus/Cow/Chicken/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	3921
Gene Symbol	RPSA
UniProt ID	RSSA_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---