

## Anti-TAp63 antibody [TAp63-4.1] (STJ16103399)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IP
Host / Source	Mouse
Reactivity	Mouse/Rat/Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS) , pH 7.4, 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---