

## Anti-Tenascin C antibody [T2H5] (STJ16103400)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P/IP/WB        |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                         |
|---------------------|---------------------------------------------------------|
| Clonality           | Monoclonal                                              |
| Concentration       | 1 mg/mL                                                 |
| Conjugation         | Unconjugated                                            |
| Purification        | Purified by protein-A affinity chromatography.          |
| Formulation         | Tris buffered saline (TBS) , pH 8.0, 15 mM sodium azide |
| Isotype             | IgG1                                                    |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.          |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|