

Anti-Vimentin antibody [VI-10] (STJ16103408)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IP/WB/IHC-P/ICC
Host / Source	Mouse
Reactivity	Pig/Chicken/Rat/Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Formulation	Tris buffered saline (TBS) , pH 8.0, 15 mM sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	7431
Gene Symbol	VIM
UniProt ID	VIME_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---