

Mouse Anti-Human IgM HRP antibody {HRP} (STJ16103413)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	HRP
Purification	Purified antibody is conjugated with activated horseradish peroxidase (HRP) of high specific enzyme activity under optimum conditions and unconjugated antibody and free HRP are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS) , pH 7.4, 0.01% thimerosal
Isotype	IgG5
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---