

## Anti-Recombinant-KSP-Cadherin, recombinant antibody (STJ16103414)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	Unconjugated
Purification	Purified protein A/G
Dilution Range	IHC-P 1:100-1:200
Formulation	10mM PBS with 0.05% BSA & 0.05% azide
Isotype	IgG
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---