

## Anti-FIBRINOGEN antibody (STJ16103416)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	1.0mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	Phosphate buffered saline<0.1% Sodium Azide (NaN <sub>3</sub> )
Storage Instruction	Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---