

## Anti-GAA antibody [2G7] (STJ16103428)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/IHC/ICC/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Adding 0.2 mL of distilled water =f 500 Mu g/mL.
Conjugation	Unconjugated
Purification	Lyophilized
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg Na <sub>3</sub> NA
Isotype	IgG2b
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2548</a>
Gene Symbol	<a href="#">GAA</a>
UniProt ID	<a href="#">LYAG_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------