

Anti-CD46 antibody [9E9] (STJ16103429)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Adding 0.2 mL of distilled water =f 500 Mu g/mL.
Conjugation	Unconjugated
Purification	Lyophilized
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3. NA
Isotype	IgG1
Storage Instruction	Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4179
Gene Symbol	CD46
UniProt ID	MCP_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---