

**Anti-PP2A-alpha/PPP2CA antibody [3B6] (STJ16103449)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC/ICC/IF/FC
Host / Source	Mouse
Reactivity	Human/Monkey/Mouse/Rat

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	Adding 0.2 mL of distilled water =f 500 Mu g/mL.
Conjugation	Unconjugated
Purification	Lyophilized
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg NaN <sub>3</sub> . NA
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---