

## Anti-MCM2 antibody [11C4] (STJ16103451)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC/IF/FC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Adding 0.2 mL of distilled water =f 500 Mu g/mL.
Conjugation	Unconjugated
Purification	Lyophilized
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3. NA
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">4171</a>
Gene Symbol	<a href="#">MCM2</a>
UniProt ID	<a href="#">MCM2_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---