

## Anti-PCK2 antibody [3F7] (STJ16103462)

### GENERAL INFORMATION

|               |                        |
|---------------|------------------------|
| Product Type  | Primary antibodies     |
| Applications  | WB/IHC/ICC/IF/FC       |
| Host / Source | Mouse                  |
| Reactivity    | Human/Mouse/Rat/Monkey |

### PRODUCT PROPERTIES

|                     |                                                                              |
|---------------------|------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                   |
| Concentration       | Adding 0.2 mL of distilled water =f 500 Mu g/mL.                             |
| Conjugation         | Unconjugated                                                                 |
| Purification        | Lyophilized                                                                  |
| Formulation         | Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3. NA |
| Isotype             | IgG2b                                                                        |
| Storage Instruction | Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.                 |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">5106</a>        |
| Gene Symbol | <a href="#">PCK2</a>        |
| UniProt ID  | <a href="#">PCKGM_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|