

**Anti-PARP/PARP1 antibody [10E11] (STJ16103477)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC/ICC/IF/FC   |
| Host / Source | Mouse              |
| Reactivity    | Human/Mouse/Rat    |

**PRODUCT PROPERTIES**

|                     |                                                                              |
|---------------------|------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                   |
| Concentration       | Adding 0.2 mL of distilled water =f 500 Mu g/mL.                             |
| Conjugation         | Unconjugated                                                                 |
| Purification        | Lyophilized                                                                  |
| Formulation         | Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3. NA |
| Isotype             | IgG2a                                                                        |
| Storage Instruction | Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.                 |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|