

## Anti-MCM5 antibody [4G10] (STJ16103486)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC/IF/FC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Adding 0.2 mL of distilled water =f 500 Mu g/mL.
Conjugation	Unconjugated
Purification	Lyophilized
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg NaN <sub>3</sub> . NA
Isotype	IgG2a
Storage Instruction	Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">4174</a>
Gene Symbol	<a href="#">MCM5</a>
UniProt ID	<a href="#">MCM5_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---