

## Anti-SAMHD1 antibody [10H8] (STJ16103495)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC/ICC/IF/FC   |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                              |
|---------------------|------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                   |
| Concentration       | Adding 0.2 mL of distilled water =f 500 Mu g/mL.                             |
| Conjugation         | Unconjugated                                                                 |
| Purification        | Lyophilized                                                                  |
| Formulation         | Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3. NA |
| Isotype             | IgG2b                                                                        |
| Storage Instruction | Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.                 |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">25939</a>       |
| Gene Symbol | <a href="#">SAMHD1</a>      |
| UniProt ID  | <a href="#">SAMH1_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|