

Anti-IGF2R antibody [6G2] (STJ16103511)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Adding 0.2 mL of distilled water =f 500 Mu g/mL.
Conjugation	Unconjugated
Purification	Lyophilized
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg Na ₃ NA
Isotype	IgG2b
Storage Instruction	Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3482
Gene Symbol	IGF2R
UniProt ID	MPRI_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---