

Anti-PGP9.5 antibody [3E4] (STJ16103512)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/FC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Adding 0.2 mL of distilled water =f 500 Mu g/mL.
Conjugation	Unconjugated
Purification	Lyophilized
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ . NA
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7345
Gene Symbol	UCHL1
UniProt ID	UCHL1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---