

Anti-Rab11A antibody [4H9] (STJ16103522)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | WB/IHC/ICC/IF/FC |
| Host / Source | Mouse |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | Adding 0.2 mL of distilled water =f 500 Mu g/mL. |
| Conjugation | Unconjugated |
| Purification | Lyophilized |
| Formulation | Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ . NA |
| Isotype | IgG2b |
| Storage Instruction | Store at -20°C for up to 1 year. Minimise freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 8766 |
| Gene Symbol | RAB11A |
| UniProt ID | RB11A_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|