

Anti-splicing factor 1 antibody [2F5D10] (STJ16103542)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC/IF/FC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Adding 0.2 mL of distilled water =f 500 Mu g/mL.
Conjugation	Unconjugated
Purification	Lyophilized
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3. NA
Isotype	IgG2a
Storage Instruction	Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7536
Gene Symbol	SF1
UniProt ID	SF01_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---