

Anti-YWHAE antibody [3G11G2] (STJ16103572)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/IHC/ICC/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Adding 0.2 mL of distilled water =f 500 Mu g/mL.
Conjugation	Unconjugated
Purification	Lyophilized
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3. NA
Isotype	IgG2b
Storage Instruction	Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7531
Gene Symbol	YWHAE
UniProt ID	1433E_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---