

Anti-KLHL12 antibody [2G11D1] (STJ16103627)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Adding 0.2 mL of distilled water =f 500 Mu g/mL.
Conjugation	Unconjugated
Purification	Lyophilized
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3. NA
Isotype	IgG2a
Storage Instruction	Store at-20°C for up to 1 year. Minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	59349
Gene Symbol	KLHL12
UniProt ID	KLH12_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---