

Anti-Keratin 8 antibody [C-51] (STJ16103633)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/IP/WB/ICC
Host / Source	Mouse
Reactivity	Pig/Cow/Human/Sheep

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---