

## Anti-Keratin 18 antibody [C-04] (STJ16103634)

### GENERAL INFORMATION

|               |                          |
|---------------|--------------------------|
| Product Type  | Primary antibodies       |
| Applications  | IHC-P/FC/IP/WB/ICC/ELISA |
| Host / Source | Mouse                    |
| Reactivity    | Mammalian                |

### PRODUCT PROPERTIES

|                     |                                                                  |
|---------------------|------------------------------------------------------------------|
| Clonality           | Monoclonal                                                       |
| Concentration       | 1 mg/mL                                                          |
| Conjugation         | Unconjugated                                                     |
| Formulation         | Phosphate buffered saline (PBS) solution with 15 mM sodium azide |
| Isotype             | IgG1                                                             |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.                   |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|