

Anti-Keratins 4, 5, 6, 8, 10, 13 and 18 antibody [C-11] (STJ16103637)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/FC/IP/IHC-P/ICC
Host / Source	Mouse
Reactivity	Mammalian

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---