

Anti-Keratin 9 antibody [Ks 9.70/Ks 9.216] (STJ16103644)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ICC/IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Supernatant
Dilution Range	ICC 1:10 IF 1:10 IHC-F 1:20 IHC-P 1:20 WB 1:25-1:300
Formulation	Hybridoma cell culture supernatant containing 0.09% sodium azide
Isotype	IgG1/IgG3
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---