

Anti-IgG1 antibody [LO-MG1-2] (STJ16103663)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/ELISA
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.1% sodium azide
Isotype	IgGk
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

TARGET INFORMATION

Gene Symbol	NA Uniprot ID= IGG1_HUMAN "
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---