

## Anti-IgG3 antibody [LO-MG3-7] (STJ16103670)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/ELISA
Host / Source	Rat
Reactivity	Mouse/Swine/Horse/Dog

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.1% sodium azide
Isotype	IgMk
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---