

Anti-IgM antibody [LO-MM-8] (STJ16103677)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/ELISA
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.1% sodium azide
Isotype	IgG2ak
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

TARGET INFORMATION

Gene Symbol	NA Uniprot ID= IGM_HUMAN "
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---