

Anti-IgG2a antibody [LO-MG-2a9] (STJ16103679)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-F/ELISA
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	PBS pH 7.2/50% (v/v) Glycerol
Isotype	IgG1k
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---