

## Anti-Mouse-IgG2a Kappa antibody [LO-MK-1] (STJ16103682)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ELISA/FC/IF/WB
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.1% sodium azide
Isotype	IgG2ak
Storage Instruction	Store at -20°C for up to 1 year. Minimise freeze-thaw cycles.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---